



MEDI-FLO®
ELASTOMERIC PUMP

Product Guide

Disposable Elastomeric
Infusion Pumps for Short-Term
and Long-Term Infusions

Delivering Efficiency to Healthcare®

Distributed by: NDC Homecare, LLC
415 New Sanford Road / La Vergne, TN 37086
www.preferredmedical.com / 615.722.7865



MEDI-FLO®

ELASTOMERIC PUMP

Why Wolf-Pak® MEDI-FLO® Elastomeric Pumps?

MEDI-FLO® Elastomeric Pumps are proven to be safe and effective and are a convenient alternative to electronic infusion pumps.

Ideally suited for home infusion, long-term care, and outpatient chemotherapy treatments. MEDI-FLO® pumps give the patient mobility and freedom to maintain an active lifestyle.

As the pharmacist's preferred disposable pump, MEDI-FLO® is easy to fill and color-coded for quick and accurate device identification. A wide variety of sizes and flow rates are available, along with an extensive library of drug stability data.



No batteries, power cords or IV poles required



Fixed flow rate design requires no programming



Guaranteed savings and easy to use



Widest product range, configurations and offering



1

Tethered cap avoids misplacing luer cap after filling device.

2

Color-coded fill port for easy identification of flow rate.

3

Dual-layered drug reservoir membrane provides positive pressure for accurate flow rate.

4

Soft material provides ultimate patient comfort.



5

Flow restrictor reduces temperature effects and maintains flow rate accuracy with skin contact.



7

GVS 1.2 micron filter provides superior performance and contamination protection. Color-coded for flow rate identification.



6

Ergonomical clamp is easy to close and open with one hand.



8

Short and long duration pumps for various therapies such as antibiotics, chemotherapy and pain management.



	Item Code	Nominal Fill	Nominal Flow Rate	Max Fill	Duration	Qty/Case	Color Code
SHORT DURATION	MDE050500	50 mL	50 mL/hr	60 mL	1 hr	48/cs	Dark Green
	MDE100500	100 mL	50 mL/hr	110 mL	1 hr	48/cs	Dark Green
	MDE101000	100 mL	100 mL/hr	120 mL	1 hr	48/cs	White
	MDE102000	100 mL	200 mL/hr	100 mL	0.5 hr	48/cs	Light Blue
	MDE200100	200 mL	100 mL/hr	300 mL	2 hr	24/cs	White
	MDE200200	200 mL	200 mL/hr	300 mL	1 hr	24/cs	Light Blue
	MDE250050	250 mL	50 mL/hr	275 mL	5 hr	24/cs	Dark Green
	MDE250100	250 mL	100 mL/hr	300 mL	2.5 hr	24/cs	White
	MDE250125	250 mL	125 mL/hr	300 mL	2 hr	24/cs	Dark Yellow
	MDE250500	250 mL	500 mL/hr	300 mL	0.5 hr	24/cs	Pink
	MDE251750	250 mL	175 mL/hr	300 mL	1.5 hr	24/cs	Dark Blue
	MDE252500	250 mL	250 mL/hr	300 mL	1 hr	24/cs	Gray
	MDE271750	270 mL	175 mL/hr	340 mL	1.5 hr	24/cs	Dark Blue
	MDE401000	400 mL	100 mL/hr	550 mL	4 hr	24/cs	White
	MDE402000	400 mL	200 mL/hr	500 mL	2 hr	24/cs	Light Blue
	MDE502500	500 mL	250 mL/hr	550 mL	2 hr	24/cs	Gray
LONG DURATION	MDL100002	100 mL	2 mL/hr	120 mL	50 hr	24/cs	Light Yellow
	MDL270002	270 mL	2 mL/hr	270 mL	135 hr	24/cs	Light Yellow
	MDL270005	270 mL	5 mL/hr	360 mL	54 hr	24/cs	Burgundy
	MDL270010	270 mL	10 mL/hr	360 mL	27 hr	24/cs	Light Green
	MDL500020	500 mL	20 mL/hr	500 mL	25 hr	24/cs	Dark Purple
CHEMOTHERAPY	MDC100020	100 mL	2 mL/hr	120 mL	50 hr	24/cs	Light Yellow
	MDC270020	270 mL	2 mL/hr	300 mL	135 hr	24/cs	Light Yellow
	MDC270050	270 mL	5 mL/hr	360 mL	54 hr	24/cs	Burgundy
	MDC270010	270 mL	10 mL/hr	360 mL	27 hr	24/cs	Light Green
	MDC300060	300 mL	6 mL/hr	360 mL	62 hr	24/cs	Dark Pink

Accessories

ECP0100	Elastomeric Carry Pouch, 5-100mL, 1/ea
ECP0250	Elastomeric Carry Pouch, 250mL, 1/ea
MDH01	MEDI-FLO® Manual Pump Filler



IN STOCK & READY TO SHIP
615-722-7865 | www.wolf-pak.com



Delivering Efficiency to Healthcare®

*Distributed By:
NDC Homecare, LLC
415 New Sanford Road
La Vergne, TN 37086
1-800-722-7865*

REV04

© 2022 Preferred Medical
www.wolf-pak.com