


Why Wolf-Pak® Biosite Plus® Needleless Valves?

Our Wolf-Pak® Biosite Plus® needleless valves feature neutral displacement connector technology, which reduces the risk of bacterial contamination and promotes patient safety.

- A single protocol is to be used with all patients, making it a seamless transition to purchase and implement clinically
- For use on all vascular catheters for blood draws or administration of IV medications
- A mechanically and microbiologically closed system reduces the risk of infection
- MRI compatible
- Improves hospital efficiencies and helps lower IV therapy costs
- Promotes aseptic technique
- Contains no latex and is non-DEHP
- The opaque valve is pressure rated for use on PICC lines
- The clear valve housing allows for visualizing the internal fluid path when flushing the connector



Item	Description	Quantity
BI8000	Clear, 60 psi	100/bx
BI8500	Opaque, 300 psi	100/bx

**CONTACT US TODAY TO EXPERIENCE
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | www.preferredmedical.com

Delivering Efficiency to Healthcare®

402 BNA Drive, Suite 500 / Nashville, TN 37217
www.preferredmedical.com / 615.722.7865