


## Why Wolf-Pak® Hypodermoclysis Kits?

Wolf Medical Supply has created easy-to-use kits to administer Hypodermoclysis Hydration Therapy. All kits contain a Wolf-Pak® Rate Flow IV Tubing, a Subcutaneous Safety Needle Set and a procedure prep kit. The dial on the rate flow tubing can be adjusted from 0 – 250 mL/hr and the set is 92" long. High Flow and Low Flow options available with either 24G or 26G bifurcated Sub-Q Needle Sets that utilize enhanced needle design and a custom luer locking mechanism as well as other patient preferred benefits.

## Why did we choose these Sub-Q Needle Sets for our Hypodermoclysis Kits?

- Sub-Q needles were designed with a patient's safety and comfort in mind
- Backcut needle tip
  - Minimizes tissue damage, scarring and local site reactions
  - Promotes faster healing
- Custom luer provides even flow to every site
- Flexible flying hinge butterfly wings fit to patients body and keep needle in place without tenting



**CONTACT US TODAY TO EXPERIENCE  
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | [www.preferredmedical.com](http://www.preferredmedical.com)

*Delivering Efficiency to Healthcare®*

402 BNA Drive, Suite 500 / Nashville, TN 37217  
[www.preferredmedical.com](http://www.preferredmedical.com) / 615.722.7865



# WOLF-PAK® HYPODERMOCLYSIS KITS

## Kits Include:

- (1) Bifurcated Sub-Q Set
- (2) Tegaderm Dressings 2" x 3"
- (1) Wolf-Pak® Rate Flow Tubing Set 92",  
0 - 250mL/hr dial
- (1) Sterile Prep Kit that includes:
  - (1) Pair of Gloves Latex Free
  - (2) Alcohol Prep Pads
  - (2) Prevantics Swabs
  - (1) Roll of Tape



| Item      | Description                           | Quantity |
|-----------|---------------------------------------|----------|
| WPSQ22406 | High Flow Kit with 24G x 6 mm Needles | 5/cs     |
| WPSQ22609 | Low Flow Kit with 26G x 9 mm Needles  | 5/cs     |

**CONTACT US TODAY TO EXPERIENCE  
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | [www.preferredmedical.com](http://www.preferredmedical.com)

*Delivering Efficiency to Healthcare®*

402 BNA Drive, Suite 500 / Nashville, TN 37217  
[www.preferredmedical.com](http://www.preferredmedical.com) / 615.722.7865