



Welcome to our Wolf-Pak® Needle-Free Family

With twenty years of experience in the infusion market, Preferred Medical invites you to join our needle-free family! Our four uniquely designed, luer activated Wolf-Pak® needleless valves unite to form a complete offering of safe and trusted product.

Whether it is displacement type, color preference or simply price – we have an option (or maybe two!) just for you. Our dedicated and knowledgeable team is available to guide you to the valve that caters to your specific needs.

SmartSite® and Extension Sets

Our Wolf-Pak® SmartSite® needleless valves features passive displacement technology and is one of the most widely used needleless valves with a proven reputation as a reliable system. The Wolf-Pak® SmartSite® offers a flat, smooth surface to avoid bacteria entry.



Item	Description	Quantity
SV2000	Needleless Valve	100/bx
SV2002	6" Extension Set	100/bx
SV2205	20mm Vial Adapter	100/bx

Biosite Plus® Clear and Extension Sets

Our Wolf-Pak® Biosite needleless valves feature neutral displacement connector technology, which reduces the risk of bacterial contamination and promotes patient safety. The clear housing allows for visualizing the internal fluid path when flushing the connector.



Item	Description	Quantity
BI8000	Clear Needleless Valve	100/bx
BI8002	Clear w/8" Extension Set	100/bx
BI8004	Clear w/14" Extension Set	100/bx

Wolf-Site® and Extension Sets

Our Wolf-Pak® Wolf-Site® needleless valves features neutral displacement connector technology and has a non-restricted high flow rate due to a unique luer activated split septum design.



Item	Description	Quantity
WS2000	Needleless Valve	100/bx
WS2002	8" Extension Set	100/bx
WS2004	14" Extension Set	100/bx

Biosite Plus® Opaque and Extension Sets

Our Wolf-Pak® Biosite needleless valves feature neutral displacement connector technology, which reduces the risk of bacterial contamination and promotes patient safety. The opaque Biosite Plus® valve is pressure rated for use on PICC lines.



Item	Description	Quantity
BI8500	Opaque Needleless Valve	100/bx
BI8502	Opaque w/8" Extension Set	100/bx
BI8504	Opaque w/14" Extension Set	100/bx



	SmartSite®	Wolf-Site®	Biosite Plus® Clear	Biosite Plus® Opaque
Displacement	Passive	Neutral	Neutral	Neutral
Priming Volume (mL)	0.11	0.09	0.04	0.04
Pressure Rating	325 psi	360 psi	120 psi	350 psi
Device Color	White top half, clear bottom half, blue internal mechanism	Clear outer housing, blue internal mechanism	Clear	Opaque
Type of Connection	Luer lock with IV set or syringe	Split septum accessed with luer of IV set or syringe	Luer lock with IV set or syringe	Luer lock with IV set or syringe
Internal Mechanism	Mechanical valve	None	Mechanical valve with internal blunt cannula	Mechanical valve with internal blunt cannula
MRI Use	Yes	Yes	Yes	Yes
Power Injectable	Yes	Yes	No	Yes
Clamping	Before IV set or syringe disconnection	Before IV set or syringe disconnection	Before or after IV set or syringe disconnection	Before or after IV set or syringe disconnection

**CONTACT US TODAY TO EXPERIENCE
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | www.preferredmedical.com

Delivering Efficiency to Healthcare®

402 BNA Drive, Suite 500 / Nashville, TN 37217
www.preferredmedical.com / 615.722.7865