

Why Wolf-Pak® TPN Bags?

Wolf-Pak® TPN Bags are made of durable EVA material (ethylene vinyl acetate) which reduces concerns about PVC extractables interacting with the lipids of TPN solutions. Our TPN Bags also eliminate the need for separate sealing clips with their detachable transfer sets & pre-attached non-reopening clamp and feature a durable reinforced hanger.


Item	Description	Quantity
Gravity/Pooling Bags		
TPNGR500ML	Gravity, 500 mL, 3 leg	50/cs
TPNGR1L	Gravity, 1,000 mL, 3 leg	50/cs
TPNGR2L	Gravity, 2,000 mL, 3 leg	50/cs
TPNGR3L	Gravity, 3,000 mL, 3 leg	50/cs

Compounder Bags

TPNCP1L	Compounder, 1,000 mL	50/cs
TPNCP2L	Compounder, 2,000 mL	50/cs
TPNCP3L	Compounder, 3,000 mL	50/cs

Compounder Bags Exactamix

TPNEX1L	Compounder, 1,000 mL	50/cs
TPNEX2L	Compounder, 2,000 mL	50/cs
TPNEX3L	Compounder, 3,000 mL	50/cs
TPNEX4L	Compounder, 4,000 mL	50/cs



**CONTACT US TODAY TO EXPERIENCE
UNBEATABLE VALUE & SERVICE**

(800) 722-7865 | www.preferredmedical.com

Delivering Efficiency to Healthcare®

402 BNA Drive, Suite 500 / Nashville, TN 37217
www.preferredmedical.com / 615.722.7865