



# **MEDI-FLO<sup>®</sup>** ELASTOMERIC PUMP

## **Product Guide**

Disposable Elastomeric Infusion Pumps  
for Short-Term and Long-Term Infusions

*Delivering Efficiency to Healthcare<sup>®</sup>*

Distributed by: Wolf Medical Supply, An NDC Company  
13951 NW 8th Street / Sunrise, FL 33325  
[www.medi-flo.com](http://www.medi-flo.com) / 800.335.9653



# MEDI-FLO<sup>®</sup>

ELASTOMERIC PUMP

**IN STOCK & READY TO SHIP**  
(800) 335-9653 | [www.medi-flo.com](http://www.medi-flo.com)

## Why Wolf-Pak MEDI-FLO<sup>®</sup> Elastomeric Pumps?

MEDI-FLO<sup>®</sup> Elastomeric Pumps are proven to be safe and effective and are a convenient alternative to electronic infusion pumps.

Ideally suited for home infusion, long-term care, and outpatient chemotherapy treatments. MEDI-FLO<sup>®</sup> pumps give the patient mobility and freedom to maintain an active lifestyle.

As the pharmacist's preferred disposable pump, MEDI-FLO<sup>®</sup> is easy to fill and color-coded for quick and accurate device identification. A wide variety of sizes and flow rates are available, along with an extensive library of drug stability data.



**No batteries, power cords or IV poles required**



**Fixed flow rate design requires no programming**



**Guaranteed savings and easy to use**



**Widest product range, configurations and offering**



**1**  
**Tethered cap** avoids misplacing luer cap after filling device.

**2**  
**Color-coded fill port** for easy identification of flow rate.

**3**  
**Dual-layered drug reservoir membrane** provides positive pressure for accurate flow rate.

**4**  
**Soft material** provides ultimate patient comfort.



**5**  
**Flow restrictor** reduces temperature effects and maintains flow rate accuracy with skin contact.



**7**  
**GVS 1.2 micron filter** provides superior performance and contamination protection. Color-coded for flow rate identification.



**6**  
**Ergonomical clamp** is easy to close and open with one hand.



**8**  
**Short and long duration pumps** for various therapies such as antibiotics, chemotherapy and pain management.



|                | Item Code | Nominal Fill | Nominal Flow Rate | Max Fill | Duration | Qty/Case | Color Code   |
|----------------|-----------|--------------|-------------------|----------|----------|----------|--------------|
| SHORT DURATION | MDE050500 | 50 mL        | 50 mL/hr          | 60 mL    | 1 hr     | 48/cs    | Dark Green   |
|                | MDE100500 | 100 mL       | 50 mL/hr          | 110 mL   | 1 hr     | 48/cs    | Dark Green   |
|                | MDE101000 | 100 mL       | 100 mL/hr         | 120 mL   | 1 hr     | 48/cs    | White        |
|                | MDE102000 | 100 mL       | 200 mL/hr         | 100 mL   | 0.5 hr   | 48/cs    | Light Blue   |
|                | MDE200100 | 200 mL       | 100 mL/hr         | 300 mL   | 2 hr     | 24/cs    | White        |
|                | MDE200200 | 200 mL       | 200 mL/hr         | 300 mL   | 1 hr     | 24/cs    | Light Blue   |
|                | MDE250050 | 250 mL       | 50 mL/hr          | 275 mL   | 5 hr     | 24/cs    | Dark Green   |
|                | MDE250100 | 250 mL       | 100 mL/hr         | 300 mL   | 2.5 hr   | 24/cs    | White        |
|                | MDE250125 | 250 mL       | 125 mL/hr         | 300 mL   | 2 hr     | 24/cs    | Dark Yellow  |
|                | MDE250500 | 250 mL       | 500 mL/hr         | 300 mL   | 0.5 hr   | 24/cs    | Pink         |
|                | MDE251750 | 250 mL       | 175 mL/hr         | 300 mL   | 1.5 hr   | 24/cs    | Dark Blue    |
|                | MDE252500 | 250 mL       | 250 mL/hr         | 300 mL   | 1 hr     | 24/cs    | Gray         |
|                | MDE271750 | 270 mL       | 175 mL/hr         | 340 mL   | 1.5 hr   | 24/cs    | Dark Blue    |
|                | MDE401000 | 400 mL       | 100 mL/hr         | 550 mL   | 4 hr     | 24/cs    | White        |
|                | MDE402000 | 400 mL       | 200 mL/hr         | 500 mL   | 2 hr     | 24/cs    | Light Blue   |
| MDE502500      | 500 mL    | 250 mL/hr    | 550 mL            | 2 hr     | 24/cs    | Gray     |              |
| LONG DURATION  | MDL100002 | 100 mL       | 2 mL/hr           | 120 mL   | 50 hr    | 24/cs    | Light Yellow |
|                | MDL270002 | 270 mL       | 2 mL/hr           | 270 mL   | 135 hr   | 24/cs    | Light Yellow |
|                | MDL270005 | 270 mL       | 5 mL/hr           | 360 mL   | 54 hr    | 24/cs    | Burgundy     |
|                | MDL270010 | 270 mL       | 10 mL/hr          | 360 mL   | 27 hr    | 24/cs    | Light Green  |
|                | MDL500020 | 500 mL       | 20 mL/hr          | 500 mL   | 25 hr    | 24/cs    | Dark Purple  |
| CHEMOTHERAPY   | MDC100020 | 100 mL       | 2 mL/hr           | 120 mL   | 50 hr    | 24/cs    | Light Yellow |
|                | MDC270020 | 270 mL       | 2 mL/hr           | 300 mL   | 135 hr   | 24/cs    | Light Yellow |
|                | MDC270050 | 270 mL       | 5 mL/hr           | 360 mL   | 54 hr    | 24/cs    | Burgundy     |
|                | MDC270010 | 270 mL       | 10 mL/hr          | 360 mL   | 27 hr    | 24/cs    | Light Green  |
|                | MDC300060 | 300 mL       | 6 mL/hr           | 360 mL   | 62 hr    | 24/cs    | Dark Pink    |

### Accessories

|         |  |
|---------|--|
| ECP0100 | Elastomeric Carry Pouch, 5-100mL, 1/ea |
| ECP0250 | Elastomeric Carry Pouch, 250mL, 1/ea   |
| MDH01   | MEDI-FLO® Manual Pump Filler           |



**IN STOCK & READY TO SHIP**  
 (800) 335-9653 | [www.medi-flo.com](http://www.medi-flo.com)

# *Delivering Efficiency to Healthcare<sup>®</sup>*

*Distributed by:  
Wolf Medical Supply, An NDC Company  
13951 NW 8th Street  
Sunrise, FL 33325  
800.335.9653*

REV03

© 2020 Wolf Medical Supply, An NDC Company  
[www.wolfmed.com](http://www.wolfmed.com)