



# WOLF-PAK® TPN BAGS

## Why Wolf-Pak® TPN Bags?

Wolf-Pak® TPN Bags are made of durable EVA material (ethylene vinyl acetate) which reduces concerns about PVC extractables interacting with the lipids of TPN solutions. Our TPN Bags also eliminate the need for separate sealing clips with their detachable transfer sets & pre-attached non-reopening clamp and feature a durable reinforced hanger.

Item	Description	Quantity
<i>Pooling Bags</i>		
TPNGR500ML	Gravity, 500 mL, 3 leg	50/cs
TPNGR1L	Gravity, 1,000 mL, 3 leg	50/cs
TPNGR2L	Gravity, 2,000 mL, 3 leg	50/cs
TPNGR3L	Gravity, 3,000 mL, 3 leg	50/cs

### Compounding Machine Compatible

TPNCP1L	Compounder, 1,000 mL	50/cs
TPNCP2L	Compounder, 2,000 mL	50/cs
TPNCP3L	Compounder, 3,000 mL	50/cs

### Exacta Mix Compatible Only

TPNEX1L	Compounder, 1,000 mL	50/cs
TPNEX2L	Compounder, 2,000 mL	50/cs
TPNEX3L	Compounder, 3,000 mL	50/cs
TPNEX4L	Compounder, 4,000 mL	50/cs

**IN STOCK & READY TO SHIP**

(800) 335-9653 | [www.wolfmed.com](http://www.wolfmed.com)

*Delivering Efficiency to Healthcare®*

Distributed by: Wolf Medical Supply, An NDC Company  
13951 NW 8th Street | Sunrise, FL 33325  
[www.wolfmed.com](http://www.wolfmed.com) | 800.335.9653



TPNGR500ML



TPNCP1L

TPNEX4L