




WOLF-PAK® HYPODERMOCLYSIS KIT

Why Wolf-Pak® Hypodermoclysis Kit?

Wolf Medical Supply has created an easy-to-use kit to administer Hypodermoclysis Hydration Therapy. The kit contains a Wolf-Pak® Rate Flow IV Tubing, a Subcutaneous Safety Needle Set, and a procedure prep kit. The dial on the rate flow tubing can be adjusted from 0 - 250mL/hr and the set is 92" long. The Sub-Q Needle Set is a 24 gauge, bifurcated set that utilizes an enhanced needle design and custom luer locking mechanism as well as other patient preferred benefits.

Why did we choose this Sub-Q Needle Set for our Hypodermoclysis Kit?

- Sub-Q needles were designed with a patient's safety and comfort in mind.
- Backcut needle tip.
 - Minimizes tissue damage, scarring and local site reactions.
 - Promotes faster healing.
- Custom luer provides even flow to every site.
- Flexible flying hinge butterfly wings fit to patients body and keep needle in place without tenting.



IN STOCK & READY TO SHIP

(800) 335-9653 | www.wolfmed.com

Delivering Efficiency to Healthcare®

Distributed by: Wolf Medical Supply, An NDC Company
13951 NW 8th Street / Sunrise, FL 33325
www.wolfmed.com / 800.335.9653

REV03



WOLF-PAK® HYPODERMOCLYSIS KIT

WPSQ22406 Kit Includes:

- (1) Bifurcated Sub-Q Set 24G x 6mm
- (2) Tegaderm Dressings 2" x 3"
- (1) Wolf-Pak® Rate Flow Tubing Set 92", 0 - 250mL/hr dial
- (1) Sterile Prep Kit that includes:
 - (1) Pair of Gloves Latex Free
 - (2) Alcohol Prep Pads
 - (2) Preventics Swabs
 - (1) Roll of Tape



Item	Description	Quantity
WPSQ22406	Hypodermoclysis Kit	5/cs

IN STOCK & READY TO SHIP

(800) 335-9653 | www.wolfmed.com

Delivering Efficiency to Healthcare®

Distributed by: Wolf Medical Supply, An NDC Company
13951 NW 8th Street / Sunrise, FL 33325
www.wolfmed.com / 800.335.9653

REV03